



CryoSuccess

The CryoSuccess is a Swiss made cryosurgery device used to liquid freeze with dexterity. The device allows direct contact between the borosilicate glass tip and the skin, allowing a constant penetrating temperature of -89°C .



N₂O Canister

300 seconds of use per canister



CryoSuccess Pen

with N₂O canister and 1mm standard tip

CryoSurgery

The CryoSuccess can be used to treat skin lesions on various animals including dogs, cats, rabbits and horses. The device comes with a 2 year warranty as standard. Additional 23g canisters of N₂O are also available for purchase.

