

# Fluid Warmer

## DOUBLE LINE IV

CODE: TTM09

**This fluid warmer uses conductive thermal transfer to warm IV fluids as it passes through the fluid administration line.**



- Compact and light weight this fluid warmer has a colourful LCD screen with detailed info
- Double groove allows for simultaneous heating of 2 IV fluid lines, and fast heat-up time means only 5minutes to reach set temperature (1)
- With a high and low temperature alarm, the unit automatically shuts off when overheating is detected (2)
- The fluid warmer can be used on a drip pole using the in-built pole clamp. The clamp can be stored when not in use to save space (3)



### Specification

|                    |                                             |
|--------------------|---------------------------------------------|
| Tube size          | 3~4.5 mm                                    |
| Heating range      | 30°C ~ 50°C, with increment of 1°C or 0.1°C |
| Accuracy           | ≤±1 °C                                      |
| Clasiffication     | Class I / type B                            |
| Power              | AC 100-240V, 50Hz/60Hz                      |
| Net weight         | 550g                                        |
| Packing dimension: | 22.0*19.0*12.1 cm                           |