

Video Otoscope



Product Parameters⁺

Model Number	VOT-01
Display	3inch LCD
Resolution	640*480(RGB)
Power	<2W
Weight	280g
Illuminance	Natural white(6leds)
Resolution	1280*720px
Image/Video function	Yes
Output	USB output
Memory type	Micro SD card 32GB
Battery type	Rechargeable lithium battery
Capacity	1350mAh
Cycle life of battery	>300 times
Battery working time	>240mins
Charging time	<120mins



Number	Length of Scope (mm)	Working Length (mm)	Max OD (mm)	Field of View	Direction of View (DOV)	Application Scenarios
BVI-32-S1	45	30	4.2	120°	0°	Ear
BVI-32-S2	110	80	4.2	120°	0°	Nose
BVI-32-S3	190	165	4.2	120°	0°	ENT
BVI-32-S4	200	175	7	120°	0°	ENT
BVI-32-S5	200	175	5.5	120°	0°	ENT
BVI-32-S6	200	175	5.5	120°	30°	ENT
BVI-32-S7	200	175	5.5	120°	70°	ENT
BVI-32-S8	200	175	6	120°	90°	ENT
BVI-32-S9	200	175	4.2	120°	0°	ENT
BVI-32-S10	190	165	3.6	120°	0°	ENT

Endoscopy System



Hysteroscope



Bronchoscope



Ureteroscope



Choledoscope



Cystoscope



Rhinolaryngoscope



Easy to use
easy to carry



Safe and reliable
non-invasive



High-definition
vision, accurate
diagnosis

