



Small Footprint X-Ray Table

Model: CXRT-1572 / Product Code: MT001

- **Tabletop length: 150cm standard** (other sizes available on request)
- **4 ways with mechanical brakes**
- **Max weight capacity for patient supportable 100 kg**
- **Table absorption thickness 5mm**
0,8mm



How CenQuip's CXRT-1572 can Enhance your Veterinary Practice

If you're looking for a reliable, versatile, and compact x-ray table for your veterinary practice, you might want to check out CenQuip's CXRT-1572. This table is designed to meet the needs of small animal veterinarians, offering high-quality imaging and easy operation.

Here are some of the key features and benefits of the CXRT-1572:

- **Small footprint:** The CXRT-1572 has a tabletop length of 150 cm, which can fit in most exam rooms. It also has a 4-way patient support table with mechanical brakes, allowing you to position your patients with ease and accuracy.
- **Four-way** floating tabletop with mechanical braking
- **Flat panel detector:** The CXRT-1572 is compatible with various flat panel detectors up to size 17x17 inches (43 x 43 cm's).
- **Removable Grid:** The CXRT-1572 comes with a provision of adding a removable grid, this makes it easier to remove the grid easily by the handle for radiographs that do not require the use of an x-ray grid.
- **Adjustable feet:** The complete table can be levelled using the adjustable feet in case the floor level of the installation area is uneven.

With the CXRT-1572, you can perform a wide range of radiographic exams for your small animal patients, from routine check-ups to emergency cases. The table is durable, easy to use, and cost-effective, making it a great investment for your practice.

To learn more about the CXRT-1572 and other CenQuip products, visit www.cenquip.com.au or contact the sales team on (02) 9679 5733.